


# Denaturing and Diluting Libraries for the NextSeq™ 500

FOR RESEARCH USE ONLY

Introduction	3
Prepare Reagents	4
Denature and Dilute Libraries	5
Denature and Dilute PhiX Control	7
Combine Library and PhiX Control	8
Next Steps	9
Technical Assistance	



This document and its contents are proprietary to Illumina, Inc. and its affiliates ("Illumina"), and are intended solely for the contractual use of its customer in connection with the use of the product(s) described herein and for no other purpose. This document and its contents shall not be used or distributed for any other purpose and/or otherwise communicated, disclosed, or reproduced in any way whatsoever without the prior written consent of Illumina. Illumina does not convey any license under its patent, trademark, copyright, or common-law rights nor similar rights of any third parties by this document.

The instructions in this document must be strictly and explicitly followed by qualified and properly trained personnel in order to ensure the proper and safe use of the product(s) described herein. All of the contents of this document must be fully read and understood prior to using such product(s).

FAILURE TO COMPLETELY READ AND EXPLICITLY FOLLOW ALL OF THE INSTRUCTIONS CONTAINED HEREIN MAY RESULT IN DAMAGE TO THE PRODUCT(S), INJURY TO PERSONS, INCLUDING TO USERS OR OTHERS, AND DAMAGE TO OTHER PROPERTY.

ILLUMINA DOES NOT ASSUME ANY LIABILITY ARISING OUT OF THE IMPROPER USE OF THE PRODUCT(S) DESCRIBED HEREIN (INCLUDING PARTS THEREOF OR SOFTWARE) OR ANY USE OF SUCH PRODUCT(S) OUTSIDE THE SCOPE OF THE EXPRESS WRITTEN LICENSES OR PERMISSIONS GRANTED BY ILLUMINA IN CONNECTION WITH CUSTOMER'S ACQUISITION OF SUCH PRODUCT(S).

**FOR RESEARCH USE ONLY**

© 2014 Illumina, Inc. All rights reserved.

**Illumina, IlluminaDx, BaseSpace, BeadArray, BeadXpress, cBot, CSpPro, DASL, DesignStudio, Eco, GAllx, Genetic Energy, Genome Analyzer, GenomeStudio, GoldenGate, HiScan, HiSeq, Infinium, iScan, iSelect, MiSeq, MiSeqDx, Nextera, NextSeq, NuPCR, SeqMonitor, Solexa, TruGenome, TruSeq, TruSight, Understand Your Genome, UYG, VeraCode**, the pumpkin orange color, and the Genetic Energy streaming bases design are trademarks of Illumina, Inc. in the U.S. and/or other countries. All other names, logos, and other trademarks are the property of their respective owners.

## Introduction

This guide explains how to denature and dilute libraries after sample preparation to prepare them for sequencing on the Illumina® NextSeq™ 500.

This guide also explains how to prepare a PhiX control and combine libraries with the PhiX control before loading them onto the reagent cartridge.



### NOTE

*The denaturation and dilution process is not necessary for all library types.*

Some library preparation protocols result in a ready-to-use normalized concentration of pooled libraries. For more information, see the documentation for the Illumina kit you used to prepare your libraries.

## Required Consumables

The following consumables are required to denature and dilute libraries.

Illumina-Supplied Consumable	Kit Name
HT1 (Hybridization Buffer)	Component of the NextSeq 500 Kit

User-Supplied Consumable	Supplier
1 N NaOH, molecular biology-grade	General lab supplier
200 mM Tris-HCl, pH 7.0	General lab supplier

The following consumables are required to denature and dilute a PhiX control.

Illumina-Supplied Consumable	Kit Name
PhiX, 10 nM	Illumina, catalog # FC-110-3002
RSB (Resuspension Buffer)	
HT1 (Hybridization Buffer)	Component of the NextSeq 500 Kit

User-Supplied Consumable	Supplier
1 N NaOH, molecular biology-grade	General lab supplier
200 mM Tris-HCl, pH 7.0	General lab supplier

## Best Practices

- ▶ *Always* prepare freshly diluted NaOH for denaturing libraries for cluster generation. This step is essential to the denaturation process.
- ▶ To prevent small pipetting errors from affecting the final NaOH concentration, prepare at least 1 ml of freshly diluted 0.2 N NaOH.

## Prepare Reagents

The following reagents are required to denature and dilute libraries for sequencing on the NextSeq.

- ▶ **0.2 N NaOH**—A fresh dilution of NaOH is required to denature libraries for cluster generation.
- ▶ **HT1**—HT1 (Hybridization Buffer) is used to dilute denatured libraries and dilute the PhiX control.
- ▶ **RSB**—RSB (Resuspension Buffer) is used to dilute the PhiX control.

### Prepare a Fresh Dilution of NaOH

- 1 Prepare 1 ml of 0.2 N NaOH by combining the following volumes in a microcentrifuge tube:
  - Laboratory-grade water (800  $\mu$ l)
  - Stock 1.0 N NaOH (200  $\mu$ l)
- 2 Invert the tube several times to mix.



#### NOTE

A fresh dilution of 0.2 N NaOH is required for denaturing libraries and preparing a PhiX control. Set aside the remaining dilution for use within the next **12 hours**. Otherwise, discard the remaining dilution of 0.2 N NaOH.

### Prepare HT1

- 1 Remove the tube of HT1 from -15°C to -25°C storage and set aside at room temperature to thaw.
- 2 When thawed, store at 2°C to 8°C until you are ready to dilute denatured libraries.

### Prepare RSB

- 1 Remove the tube of RSB from -15°C to -25°C storage and set aside at room temperature to thaw.
- 2 When thawed, store at 2°C to 8°C until you are ready to dilute denatured libraries.

## Denature and Dilute Libraries

The loading concentration for sequencing on the NextSeq 500 is approximately 3 pM, depending on library preparation and quantification methods.

This procedure includes three variations for different starting library concentrations, either 1 nM, 2 nM, or 4 nM.



### NOTE

Typically, it is important that not more than 1 mM NaOH is in the final solution after diluting with HT1. However, introducing 200 mM Tris-HCl ensures that the NaOH is fully hydrolyzed in the final solution. As a result, template hybridization is not affected even if more than 1 mM NaOH is in the final solution after diluting with HT1.

## Denature Libraries

- 1 Combine the following volumes of libraries and freshly diluted 0.2 N NaOH in a microcentrifuge tube.

Starting Library Concentration	Library	0.2 N NaOH
4 nM	5 $\mu$ l	5 $\mu$ l
2 nM	10 $\mu$ l	10 $\mu$ l
1 nM	20 $\mu$ l	20 $\mu$ l
0.5 nM	120 $\mu$ l	120 $\mu$ l

- 2 Discard the remaining dilution of 0.2 N NaOH or set aside to prepare a PhiX control within the next 12 hours.
- 3 Vortex briefly to mix the library solution, and then centrifuge at  $280 \times g$  for 1 minute.
- 4 Incubate for 5 minutes at room temperature to denature libraries into single strands.
- 5 Add the following volume of 200 mM Tris-HCl, pH 7.

Starting Library Concentration	200 mM Tris-HCl, pH 7
4 nM	5 $\mu$ l
2 nM	10 $\mu$ l
1 nM	20 $\mu$ l

- 6 Vortex briefly to mix the library solution, and then centrifuge at  $280 \times g$  for 1 minute.

## Dilute Denatured Libraries to 20 pM

- 1 Add the following volume of pre-chilled HT1 to the tube of denatured libraries.

Starting Library Concentration	Pre-chilled HT1
4 nM	985 $\mu$ l
2 nM	970 $\mu$ l
1 nM	940 $\mu$ l
0.5 nM	2640 $\mu$ l

- 2 Vortex briefly to mix the library solution, and then centrifuge at  $280 \times g$  for 1 minute. The result is a 20 pM denatured library.
- 3 Place the 20 pM libraries on ice until you are ready to proceed to final dilution.

## Dilute Denatured Libraries to 3 pM

The following steps are an example of diluting to a loading concentration, which in the following example is 3 pM. However, the appropriate loading concentration can vary depending on library preparation and quantification methods.

- 1 Dilute the denatured 20 pM library solution to a 3 pM loading concentration, as follows:
  - 450  $\mu$ l Denatured library solution (450  $\mu$ l)
  - 2550  $\mu$ l Pre-chilled HT1 (2550  $\mu$ l)The total volume is 3 ml at 3 pM.
- 2 Invert several times to mix and then pulse centrifuge the solution.
- 3 Place the solution on ice until you are ready to load libraries onto the reagent cartridge.

## Denature and Dilute PhiX Control

Use the following steps to denature and dilute a PhiX control library.



### NOTE

For storage longer than two weeks, store the PhiX library at a concentration of 4–10 nM at -15° to -25°C. When libraries are stored at 20 pM, cluster numbers tend to decrease after two weeks of storage.

### Dilute PhiX Library from 10 nM to 4 nM

- 1 Thaw one tube of 10 nM PhiX stock (10 µl/tube).
- 2 In a 2 ml Eppendorf tube, combine the following volumes:
  - 10 nM PhiX (10 µl)
  - Resuspension Buffer (15 µl)The total volume is 25 µl at 4 nM.
- 3 Vortex briefly to mix the PhiX solution, and then pulse centrifuge.



### NOTE

[Optional] Store the 4 nM PhiX library at -15° to -25°C for up to three months.

### Denature PhiX Library and Dilute to 20 pM

- 1 In a 2 ml Eppendorf tube, combine the following volumes:
  - 4 nM PhiX library (5 µl)
  - 0.2 N NaOH, freshly diluted (5 µl)
- 2 Vortex briefly to mix the PhiX library solution, and then pulse centrifuge.
- 3 Incubate for 5 minutes at room temperature to denature the PhiX library into single strands.
- 4 After incubation, vortex briefly to mix the PhiX library solution.
- 5 Centrifuge at 280 × g for 1 minute.
- 6 Add 5 µl 200 mM Tris-HCl, pH7.
- 7 Vortex briefly and then centrifuge at 280 × g for 1 minute.
- 8 Discard the remaining dilution of 0.2 N NaOH.
- 9 Add 985 µl pre-chilled HT1. The total volume is 1 ml at 20 pM.
- 10 Vortex briefly to mix the solution.



### NOTE

[Optional] Store the denatured 20 pM PhiX library up to two weeks at -15° to -25°C. After two weeks, cluster numbers tend to decrease.

## Combine Library and PhiX Control

Illumina recommends a low-concentration PhiX control spike-in at 1% for most libraries.

- 1 Combine the following volumes of denatured PhiX control library and denatured library.

	Most Libraries (1%)
Denatured and diluted PhiX control, 20 pM	5 $\mu$ l
Denatured and diluted sample library, ~ 3 pM	2995 $\mu$ l

- 2 Set the combined library and PhiX control aside on ice until you are ready to load it onto the reagent cartridge.



## Next Steps

After denaturing and diluting your libraries and preparing the optional PhiX control, you are ready to load libraries onto the thawed reagent cartridge and set up the sequencing run. For complete instructions, see the *NextSeq 500 System User Guide* (part # 15046563).

Visit the NextSeq 500 support pages on the Illumina support website for access to documentation, software downloads, frequently asked questions, and online training.

## Notes

## Technical Assistance

For technical assistance, contact Illumina Technical Support.

**Table 1** Illumina General Contact Information

<b>Illumina Website</b>	www.illumina.com
<b>Email</b>	techsupport@illumina.com

**Table 2** Illumina Customer Support Telephone Numbers

<b>Region</b>	<b>Contact Number</b>	<b>Region</b>	<b>Contact Number</b>
North America	1.800.809.4566	Italy	800.874909
Austria	0800.296575	Netherlands	0800.0223859
Belgium	0800.81102	Norway	800.16836
Denmark	80882346	Spain	900.812168
Finland	0800.918363	Sweden	020790181
France	0800.911850	Switzerland	0800.563118
Germany	0800.180.8994	United Kingdom	0800.917.0041
Ireland	1.800.812949	Other countries	+44.1799.534000

### Safety Data Sheets

Safety data sheets (SDSs) are available on the Illumina website at [www.illumina.com/msds](http://www.illumina.com/msds).

### Product Documentation

Product documentation in PDF is available for download from the Illumina website. Go to [www.illumina.com/support](http://www.illumina.com/support), select a product, then click **Documentation & Literature**.



Illumina

San Diego, California 92122 U.S.A.

+1.800.809.ILMN (4566)

+1.858.202.4566 (outside North America)

[techsupport@illumina.com](mailto:techsupport@illumina.com)

[www.illumina.com](http://www.illumina.com)